

CD300LB-CD300 molecule-like family member b Gene

Reference: PTXBC028091

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD300b; CLM-7; CLM7; CMRF35-A2; IREM-3; IREM3; TREM-5; TREM5; CMRF35-like molecule 7; CD300 antigen like family member B; immune receptor expressed on myeloid cells 3; leukocyte mono-Ig-like receptor 5; triggering receptor expressed on myeloid cells 5; CD		

For research use only.

GENE DESCRIPTION	CD300 molecule-like family member b
GENE ID	124599
ORIGIN SPECIES	Human
ACCESSIONS	BC028091
VECTOR	pDONR223
NCBI SYMBOL	CD300LB
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgcagaaggtgcaagccagagctcgggcagaaactccagagtgcatctgggatctgcattggccactggttcagatcaggcggagcaggagccgggaaggcagagccatgtggctgccccctgctctgctcctctcagcctctcaggctgttttccatccaaggcccagagctctg
tgagagccccagagcaggggtccctgacggttcaatgccactataagcaaggatgggagacacataaagtgggtggtgcccagggggtgcgtgggatacatgcaagatcctcattgaaaccagaggggtcggaagcaggagagaagagtgaccgtgtgtccatcaaggacaatcaga
aagaccgcacgttcactgtgaccatggagggtcaggcagatgacgcagatgtttactggtgtgggattgaaagaaggagacctgacctgggactcaagtgaagtgatcgttgaccagaggagcgggtccacaacagcaagctcacctaccaacagcaatggcagtggtc
atcggtctccacaagaggaaaccactacatgctcctgtattgtgaaggtgccactctgtcatctgtgctcactgccatcctctgttgaagggtctcagaggggtccctgaggagcagggggaacagcctatcatgaactctccgaaccttgactaaagacatggccacttag