



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

DNASE1L1-deoxyribonuclease I-like 1 Gene

Reference: PTXBC028092

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

DNAS1L1; DNASEX; DNL1L; G4.8; XIB; deoxyribonuclease-1-like 1; DNase I, lysosomal-like; DNase I-like 1; DNase I-like, muscle-specific; deoxyribonuclease I-like 1; deoxyribonuclease 1 like 1

For research use only.

GENE DESCRIPTION	deoxyribonuclease I-like 1
GENE ID	1774
ORIGIN SPECIES	Human
ACCESSIONS	BC028092
VECTOR	pDONR223
NCBI SYMBOL	DNASE1L1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacactaccacactgcaccticctctctcatctcgtggccaattggggcccgagggcttctgcactctgcctctcaatgccacgcgctgacactggccaaagtgggccaggagcaggtgatggacacacttaattgctggatactgctctgctgatcatatggtgctgcaggaggtgggtgactcttcggcagcgcacatccctctctctctgcagaaactcactgaattgatgctctggccctacagaccactgagacacccccctgacacacacacactacatggagacagctgctgactctctatcggfcaacaaaacacagctctgagatctactgctatcaacatgacgagatgacgtcttcttcggcgagccacttggcccgagctcttctggccacagctctcccaacgctgtgtgtgtccctctgacacacactcaaggccgtatagagaaggagctgaacccgctgctgactgtgttctgcagctgtcccgactgtcccgacagcaaacgacgctgacatcgtctgctggtctcaatggctcgtctcactgaccaaaaggcgcctggacaagctggagctcggagctacggagcttccactgggtgatggccgatggggaggacacacacatgctcggcgccagacacactgacactatgacctgcctgctgctcacggggagctggggaggtctgctgcacacggcgctgcttctgacttcccacagacttccagctcaccgaggaggagggccctcaatcactgtagaccactacccctggaggtggagctgaagctgagcaggcgcacagcgtccagcctctcagctcactctgttctgtgctatcacctctgctcccctcagctgctgcctgctgctga

More info: Please contact us for any question: contact@proteogenix.fr