

C7orf16-chromosome 7 open reading frame 16 Gene

Reference: PTXBC028094

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C7orf16; GSBS; protein phosphatase 1 regulatory subunit 17; G-substrate		
For research use only.			

GENE DESCRIPTION	chromosome 7 open reading frame 16
GENE ID	10842
ORIGIN SPECIES	Human
ACCESSIONS	BC028094
VECTOR	pDONR223
NCBI SYMBOL	C7orf16
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatgtccactgagcaaatgcagccactggaagtctcagaagacagactggacaagctagaccctcgttcagccacttagatgatcttcagaccagttcattaaggactgtgatctcaaaaagaagcctagaaagggaataatgtacaggccaccctgaatgttgagtcagacaaaa
aaaaccaaggaggaagatacaccggcgtgcacatcccaccttcataccaggtgtgtttcagaacatttaataaaagatacagatgttcaagagagacatccaaaggcgaataatgatccctgtttctcataacactgacctggaacagaaaaagccaaggagaaaaagacacacctgccc
tgcacatgtccctttgcagcagggtgtgacattgtcagggatgagagacccaaagcaatcgtggaagatgacgaaaaggatggtgacaagatagctatttaa