

GLRX2-glutaredoxin 2 Gene

Reference: PTXBC028113

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CGI-133; GRX2; glutaredoxin 2; bA101E13.1 (GRX2 glutaredoxin (thioltransferase) 2); glutaredoxin (thioltransferase) 2		

For research use only.

GENE DESCRIPTION	glutaredoxin 2
GENE ID	51022
ORIGIN SPECIES	Human
ACCESSIONS	BC028113
VECTOR	pDONR223
NCBI SYMBOL	GLRX2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagagcaatacatcatctttggagaatttagcgacggcgctgtgaaccagatccaagaacaatttctgataattgtgtgtgattttcctcttactgtacaatggcaaaaagctttccatgacatgaatgttaactataaagtgtggaactggacgtgtgaa
tatggaaaccagttccaagatgctctttacaaaatgactggtgaaagaactgttccaagaataattgtcaatggtactttattggagggtgcaactgacactcataggttcacaaagaaggaaaattgctccactagttcatcagtggtatttaaaaaaaagtaagaggaaagaatttcagtga