

RNF32-ring finger protein 32 Gene

Reference: PTXBC028120

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FKSG33; HSD15; LMBR2; RING finger protein 32		
For research use only.			

GENE DESCRIPTION	ring finger protein 32
GENE ID	140545
ORIGIN SPECIES	Human
ACCESSIONS	BC028120
VECTOR	pDONR223
NCBI SYMBOL	RNF32
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttaaaaaataagggtcactcatctaagaagataactggcagtcagtgtgtttacaagatcacattttacatgatctcaacttcgaaatcttcagttgcagatcatttaagacacaaagtaacaaagagaacaaatctctaaaaagagatacaaaagcgaataatagatactg
gacttaaaaaaacacacagtgcccaactagaagactcagaaaaagaataatgtttgatcccaaccgccgccgttgactttggcagagaagtgggctcattgggctccaccacctccactgtcatcagatgaatgggagaaggtgaaacagcgctctctcgaagggaactcc
gtgcaaccatgccccatctgtaagaagaattcagcttcctcaggtgtgtttcatgctccatgtgtccacaagcatgtcttcaggctttgaaaagttcacaaataagaaaacctgtcctctctgtaaaagaaccagatcaaacccagtgatacacgatgggcccgcctgttc
agaatcaagtggtgaccagaatccaagcctactggagaggtgtgtttagaagtggtacagaaacctgaggaaaacagtaacctccacagatgccaaagttaagaaaaattcttgaaaaaaagttcacagaatcagccaccgcctctgtgtcatacacacacattgaagag
ctctttgcagaaatcagatcagtgcttgccataa