

C10orf53-chromosome 10 open reading frame 53 Gene

Reference: PTXBC028127

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UPF0728 protein C10orf53; chromosome 10 open reading frame 53		
For research use only.			

GENE DESCRIPTION	chromosome 10 open reading frame 53
GENE ID	282966
ORIGIN SPECIES	Human
ACCESSIONS	BC028127
VECTOR	pDONR223
NCBI SYMBOL	C10orf53
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcccaagaacgcagtggtcatcctgcgtatggccctacagcgcggcaggcctaccgggtggagcaccacacctccgctgcagggcctgcaagctgtgttggccatagatggacatgaggtcatcctagagaagatagaagactggaatgtgttggactcatggtgaatgaaga
agtcattctccactgcaacattaggacttggagtttggcaagctcaccaccaagtagtgacaagaggaccaccagcagctccaggctgacatttaccagcttagcagcccgtgtgggatgaaagtctctcttgcacagttcccaaaagaccaggacctgactgtattgtattggcc
cagattggatcatgtatccatttccaaaccaatcttggacctgggttggccaggcctggatcacatgctcatttctggactagagaagagggaactcagccctactga