



Reference: PTXBC028163

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RE2; G-protein coupled receptor 161; G-protein coupled receptor RE2; G protein-coupled receptor 161		

For research use only.

GENE DESCRIPTION	G protein-coupled receptor 161
GENE ID	23432
ORIGIN SPECIES	Human
ACCESSIONS	BC028163
VECTOR	pDONR223
NCBI SYMBOL	GPR161
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr