

PPBP-pro-platelet basic protein (chemokine (C-X-C motif) ligand 7) Gene

Reference: PTXBC028217

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	B-TG1; Beta-TG; CTAP-III; CTAP3; CTAPIII; CXCL7; LA-PF4; LDGF; MDGF; NAP-2; PBP; SCYB7; TC1; TC2; TGB; TGB1; THBGB; THBGB1; platelet basic protein; C-X-C motif chemokine 7; CXC chemokine ligand 7; beta-thromboglobulin; chemokine (C-X-C motif) ligand 7; co		

For research use only.

GENE DESCRIPTION	pro-platelet basic protein (chemokine (C-X-C motif) ligand 7)
GENE ID	5473
ORIGIN SPECIES	Human
ACCESSIONS	BC028217
VECTOR	pDONR223
NCBI SYMBOL	PPBP
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagcctcagactgataccacccctctgttaacagtgcgagaccactcatgccttcaggtgctgctgcttctgtcattgctgctgactgctctgcttctccaccaaaggacaaactagagaaactggcgaaaggcaagaggaaagtctagacagtgactgtatgctgaactcc
gctgcatgtgtataaagacaacctctggaattcatccaaaaacatccaaagtftggagtgatcgggaaaggaaacccattgcaaccaagtcgaagttagatgccacactgaaggatgggaggaaaatctgcctggaccagatgctccagaatcaagaaaattgtacagaaaaaattggc
agggtgatgaatcgtctgattaa