

CYBA-cytochrome b-245, alpha polypeptide Gene

Reference: PTXBC028224

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	p22-PHOX; cytochrome b-245 light chain; cytochrome b light chain; cytochrome b(558) alpha chain; cytochrome b(558) alpha-subunit; cytochrome b, alpha polypeptide; cytochrome b-245, alpha polypeptide; cytochrome b558 subunit alpha; flavocytochrome b-558 al		

For research use only.

GENE DESCRIPTION	cytochrome b-245, alpha polypeptide
GENE ID	1535
ORIGIN SPECIES	Human
ACCESSIONS	BC028224
VECTOR	pDONR223
NCBI SYMBOL	CYBA
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggggcagatcagtgaggccatgtgggccaacgagcagggcgtgctggcgtcggcctgagtgagtgacgtcagggacgggtggaggtgcagcctggaggggtgtcccaagacccagccgggacctcgggctacttacagggtggggaagtggggcgccaggcgggcccaggc
cgggcccgggtcaggccaggaggggtgcggggaacggggcgggacccctcaggccgggctggaaggaggagttctgagacccccagtaattcccttgacagacctgcggagcggcgccccctcccccaftcccttctctcgttccccgactccgcaaggagggaagtggca
gcgaggggaagaggcgggttcagccccgtggttccgggtcaccgccccgaagccccagtgggggctccggcctgggcatcgggagaagctccctgccccctggtagccagctggcctggaggtcgtctccctgggctggggtggggaatgggctcatgccccggggc
tcagccctctcatctggagttcgcgtggagggcgtctgagatttcacaaggccggaaccacccagtgacggccctgccagcccagctgcacaggcacttctcttgggtgggatggcacagccccagtcacccctctgtgcatttctcagaggttctggtgaggcagcgg
actccattggaaatggacctctcgggaccaggccaggaggtgtgtgtgcccagcacgggtga