

MT1M-metallothionein 1M Gene

Reference: PTXBC028280

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MT-1M; MT-IM; MT1; MT1K; metallothionein-1M; metallothionein 1K; metallothionein 1Y; metallothionein 1M		

For research use only.

GENE DESCRIPTION	metallothionein 1M
GENE ID	4499
ORIGIN SPECIES	Human
ACCESSIONS	BC028280
VECTOR	pDONR223
NCBI SYMBOL	MT1M
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacccaactgctcctgcaccactgggtgtctctgcgcctgcaccggctctgcaagtgcaagagtgcaaatgcacctctgcaagaagagctgctctgctgccccgtgggctgtgccaagtgtgccacggctgtgtctgcaaggagacgttgagaactgcagctgctgtgcctga