



Antibodies and Proteins

ProteoGenix SAS
Espace Européen de l'Entreprise
15 rue de La Haye
67300 Schiltigheim, France
(0)3 90 20 54 70 - Fax +33 (0)3 88
contact@proteogenix.fr
www.ProteoGenix.science

PLSCR4-phospholipid scramblase 4 Gene

Reference: PTXBC028354

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TRA1; phospholipid scramblase 4; Ca(2+)-dependent phospholipid scramblase 4; PL scramblase 4; cell growth inhibiting protein 43		

For research use only.

GENE DESCRIPTION	phospholipid scramblase 4
GENE ID	57088
ORIGIN SPECIES	Human
ACCESSIONS	BC028354
VECTOR	pDONR223
NCBI SYMBOL	PLSCR4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

ORF1 (Only for reference, product provided without STOP codon):

atgtacagctgtgtgtccaccagccctgaacacgctcagggtgaattggaaataaacacaacacacagatcaacaggctcgtatgctctctctgaataacattctcatattttaccaggagcccccaggacagctgtccctccacctctggtctaccaggaggcttgcctatgggatactac
atgtccacagacacccagctacctctccctgtaccacagcagttgtgttgatcatcctgtgcgctatgacgtctgcgaataatcagatcgtgtcccaataacatgatgtgcagggggcaactctatgtgcaactgcctctcgtgttggaattactgttcagtttggaacaataatc
ttctcagctatttttgagccttggtgaatgatgacatgttttgaaactaataatgatatgatataataaaaacactgcacacagatgttttactgttaacccgaagacagatgactttacaggatgctctatgcctacgaactgcctctgctctcctggtctactgtatgtatggcgagaaatc
atgacaattgcagagaccttcacagctcacctctgtttctgtctgtgccccctctgcacagaagaagctggagggtgtgctgtctctgtgttaccatttgcgtttgttggcgaaatttggaaacttggaaactgtgcagggtgtgtacagatcctcaaaatgggaagaagaatattgataagaagattctg
gggcatctcacaactatgctgtgtgttcagattctgttttgaggcttcaaacctccttatggcatalcctcaatcggcagatattatccggagtggaattgttttgcacacagatgcagatgctgacacatttggacattcctccactgaacctggatgtggaagtgaagagccatgattttgga
ctctcctctcaattcatcatgtatgtatgaagaatctccacacacagcttcaagaatg

More info: Please contact us for any question: contact@proteogenix.fr