



Reference: PTXBC028363

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MCAHS; MCAHS1; MDC4; PIG-N; GPI ethanolamine phosphate transferase 1; MCD4 homolog; phosphatidylinositol-glycan biosynthesis class N protein; phosphatidylinositol glycan anchor biosynthesis class N		

For research use only.

GENE DESCRIPTION	phosphatidylinositol glycan anchor biosynthesis, class N
GENE ID	23556
ORIGIN SPECIES	Human
ACCESSIONS	BC028363
VECTOR	pDONR223
NCBI SYMBOL	PIGN
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]