

C7orf62-chromosome 7 open reading frame 62 Gene

Reference: PTXBC028365

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C7orf62; chromosome 7 open reading frame 62		
For research use only.			

GENE DESCRIPTION	chromosome 7 open reading frame 62
GENE ID	219557
ORIGIN SPECIES	Human
ACCESSIONS	BC028365
VECTOR	pDONR223
NCBI SYMBOL	C7orf62
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtccttttcggtccataaccagaaggcagcaaaaggcctttgccactggaacctctctttttctccaagtcacgtacgtacatttcaagaagagaacaacgactacacctaagaatactctctgtaggatgtttctagtgcccaagatacaagcaaatgtagaagaa
aagatgttgctgactactatgaacaaatgtttcagtcagttttgaacatcacctaggagaagcagtgacaggaattgctcctatctatccacttccattctcatatcctcagtcctccagcgacactctctacaagttcttttagattatattggccatgtcaagatgaacagatatttttattc
aacaaatgaaaattatagtcatttctataacattccaatgaggtttttatgcaatggcatgtttcagtgataaaagtccagttatgtatctcgacgatgtgacacagtcacagtccttaaggaggtcatcacagattttctcacacaaactcataaactgtcactctaccttggccagactatgaaa
gtaggcactaaaggaccaggcgataacttacaccaagttgcacctgacctactctccagaacaaatcataaagtacttgtgcaaatccgaagaattcatggacccggcaacatttataaactgtataatagaccatacacatcaactctggattctgaggtgtatggcctgtccttcacg
tttctag