

SLC4A2-solute carrier family 4, anion exchanger, member 2 (erythrocyte membrane protein band 3-like 1) Gene

Reference: PTXBC028601

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AE2; BND3L; EPB3L1; HKB3; NBND3; anion exchange protein 2; anion exchanger 2 type a; anion exchanger 2 type b1; anion exchanger 2 type b2; erythrocyte membrane protein band 3-like 1; non-erythroid band 3-like protein; solute carrier family 4 (anion exchan		

For research use only.

GENE DESCRIPTION	solute carrier family 4, anion exchanger, member 2 (erythrocyte membrane protein band 3-like 1)
GENE ID	6522
ORIGIN SPECIES	Human
ACCESSIONS	BC028601
VECTOR	pDONR223
NCBI SYMBOL	SLC4A2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr

Page 1/1