

C12orf10-chromosome 12 open reading frame 10 Gene

Reference: PTXBC028904

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	Gamm1; MST024; MSTP024; MYG; MYG1; UPF0160 protein MYG1, mitochondrial; melanocyte proliferating gene 1; melanocyte related; chromosome 12 open reading frame 10		

For research use only.

GENE DESCRIPTION	chromosome 12 open reading frame 10
GENE ID	60314
ORIGIN SPECIES	Human
ACCESSIONS	BC028904
VECTOR	pDONR223
NCBI SYMBOL	C12orf10
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtgggaccctctatgacaagatgtatgagaacttggaggaggtggatctgtggacaatggatctccagtgaggcagagggggagcctcgatgacactaccctgagtgacagagttgctgacttaactctactgggaaccaccccgaccaagacactgaggca
gggttcaagcgtgcaatggatctgggttcaagaggagttctgcagagattagatttctaccaacacagctggctgccagccccggccttgggtgaagagggccctggccagcgattccaggtggaccgaatggagagattgtggaactggcgaaaggtgcatgtccctggaaggagcat
ctctaccactggaatctgggctgtccctccagtgcccatctcttggttatctacactgaccaggtggacagtgccgaatacagtggtgccaagagagccctactatccaaagccgctgccccgccagagccatggcgggtcttctgggacgaggccctggaccaggtcagtg
ggatccctggctgcatctctccatgcaagcgggtcattggcgggtcaccgaccccgagaggggtgccttgagcatggccctggccacctggccagcgctcaccctcccaaatctcctag