

DNAJB9-DnaJ (Hsp40) homolog, subfamily B, member 9 Gene

Reference: PTXBC028912

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ERdj4; MDG-1; MDG1; MST049; MSTP049; dnaJ homolog subfamily B member 9; DnaJ (Hsp40) homolog, subfamily B, member 9; ER-resident protein ERdj4; endoplasmic reticulum DNA J domain-containing protein 4; endoplasmic reticulum DnaJ homolog 4; microvascular en		

For research use only.

GENE DESCRIPTION	DnaJ (Hsp40) homolog, subfamily B, member 9
GENE ID	4189
ORIGIN SPECIES	Human
ACCESSIONS	BC028912
VECTOR	pDONR223
NCBI SYMBOL	DNAJB9
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctactccccagtcgaatttcatcttgcatttgaataacagaatttctgacctcaaaagctactatgatacttaggtgtgccaaatcggcatcagagcgccaaatcaagaaggcctttcacaagttggccatgaagtaccacctgacaaaaataagagccggatg
ctgaagcaaaatcagagagattgcagaagcatatgaacactctcagatgctaataagacgaaaagatgatacacttggacacagtccttttactagtggttaaggacaaagaggtagtggaagttctttgagcagtcatttaactcaattttgatgactatttaagactttggctttttgg
tcaaaacaaaacactggatccaagaagcgttttgaataatcttcagacacgcccagatgggtgtccagtagacaaagcatttccaagaattttcttggagggtgattattgatgacatgtttgaagatatggagaaaatgtttcttttagtgggtttgactctacaatcagcataca
gtacagactgaaaatagatttcatggtatgtagcaagcactgcaggactgtcactcaacgaagagaaatattgttactacatacactgactgttcaggacagtag