

C10orf111-chromosome 10 open reading frame 111 Gene

Reference: PTXBC029034

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C10orf111; bA455B2.4; chromosome 10 open reading frame 111		

For research use only.

GENE DESCRIPTION	chromosome 10 open reading frame 111
GENE ID	221060
ORIGIN SPECIES	Human
ACCESSIONS	BC029034
VECTOR	pDONR223
NCBI SYMBOL	C10orf111
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaatccctgcagactcccagaccgcgaaatcaagataaaaggagagatgggttaaacaatgcctatgggcaataatgcaggaatcttgagcccgaaaagagaaaggcagtaagagttgccttgagttcagcaacagctgcacagaatcccgccagtgccca
ctgtggctgtccaagcaatggagactcaggtaccatcgagtcgtgcaaaagtcgggacaaagtgaagcggccgaataacatttaaagctcaggaatctggtgatctaccctggccagaacttagaacggcaggttgcttctgacctaatgagcctcctcttccc
cgcttttccggccttacttgggccccttcttctctttacgtatctgcctcttttcaatctcctcactgtgtgtttgtatcctattttctggtatag