

ATP4B-ATPase, H⁺/K⁺ exchanging, beta polypeptide Gene

Reference: PTXBC029059

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ATP6B; potassium-transporting ATPase subunit beta; ATPase, H ⁺ /K ⁺ exchanging, beta polypeptide; ATPase, H ⁺ /K ⁺ transporting, beta polypeptide; gastric H ⁺ /K ⁺ ATPase subunit beta; gastric H ⁺ /K ⁺ ATPase beta subunit; gastric hydrogen-potassium ATPase, beta;		

For research use only.

GENE DESCRIPTION	ATPase, H ⁺ /K ⁺ exchanging, beta polypeptide
GENE ID	496
ORIGIN SPECIES	Human
ACCESSIONS	BC029059
VECTOR	pDONR223
NCBI SYMBOL	ATP4B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgctctgcaggagaagaagacgtgtggccagcgcattggagggtccagcgttactgctggaacccggacacggggcagatgctggccgcaccctgtcccggtgggtgtggatcagcctgtactacgtggccttctactgtgtgatgactgggctcttgcctgtgcctctat
gtgctgatgcagacagtgaccgtacacaccggactaccaagaccagctacgggtaccagggttaaccttaaggccggaatgtttacggggagaaaggcctggaaatgtctacaacgtctctgataacagaacctgggcagacctcacacagactctccagcctctagcaggcta
cttccagcagcccgaggagacatcaactgcactccgagcagctactcttccaggagaggttccgctcccaaccacacaaagttctcctgcaagttcacggcagatagctgcagaactgctcagcctggcgatcccaactcggctttgaagaaggaaagccatgtttattat
taaaatgaacagatgtcaagttctctcccagcaacggctcggccccagagtgactgcgcttctggaccagccccgcagctcggccagccgctgcaggtcaagtactaccctcccaacggcacctcagctgcactacttccctattacgggaagaaagccagccccacta
cagcaacccctgggtggcagcgaagctcctcaacatccccaggaacgctgaggtcgccatcgtgtgcaaggtcagtcagcagtgaccttcaacaatccccacgaccggtatgaagggaagtgaggttcaactcaagattgagaagtga