



Reference: PTXBC029121

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD212; IL-12R-BETA1; IL12RB; IMD30; interleukin-12 receptor subunit beta-1; IL-12 receptor beta component; IL-12 receptor subunit beta-1; IL-12R subunit beta-1; cluster of differentiation 212; interleukin 12 receptor, beta 1; interleukin-12 receptor beta-		

For research use only.

GENE DESCRIPTION	interleukin 12 receptor, beta 1
GENE ID	3594
ORIGIN SPECIES	Human
ACCESSIONS	BC029121
VECTOR	pDONR223
NCBI SYMBOL	IL12RB1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]