

LYG1-lysozyme G-like 1 Gene

Reference: PTXBC029126

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SALW1939; lysozyme g-like protein 1; lysozyme G-like 1; lysozyme g1		
For research use only.			

GENE DESCRIPTION	lysozyme G-like 1
GENE ID	129530
ORIGIN SPECIES	Human
ACCESSIONS	BC029126
VECTOR	pDONR223
NCBI SYMBOL	LYG1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctgatttggctgctgctggcctctgcccctgatggactgtctgaaagcagcaactgggagctatggaacatccaaagcctggacccctggagcatcttgggattggaagacgtcacggcctgaactactgtggagttcgtgcttctgaaaggctggctgaaatagac
atgccatacctctgaaatataacccatgatgcaaaccttggccaaaagtactgcatggatctgcccgtgacgctggtgtcttgtccaggaagtctccgggtgacaaaattctgtcaacatggcgcataggactagcatggcagggacctggctctcaagctccacatcctggatta
gtgagttcagggttccagacaactgaagtctgactactagaatcaagaatccagaggaggtttccaacctggacccctgaccagtacctgagaggttgactctgtgcctacagtgggggtgctgctatgtccgaagcagcaggacctgagctgtgacttctgcaatgatgtccttg
cacgagccaagtacctaagagacatggcttctaa