

DPY19L3-dpy-19-like 3 (C. elegans) Gene

Reference: PTXBC029162

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	dpy-19-like 3; dpy-19-like protein 3; protein dpy-19 homolog 3; dpy-19 like 3 (C. elegans)		

For research use only.

GENE DESCRIPTION	dpy-19-like 3 (C. elegans)
GENE ID	147991
ORIGIN SPECIES	Human
ACCESSIONS	BC029162
VECTOR	pDONR223
NCBI SYMBOL	DPY19L3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatgccatccgcaagaagagaataagagccacagaagttctgaagacttccagcccaagaagaatgtgaagttgaaaataaattgccatctgtgtaccagtagaagattatggaagattttgattgacaattggtggaaccattgcccttgattggacttcttacatctg
tctaccttgccacgttacatgaaaatgatttatggttttctaataaaggatggagtttcttgaccattgtatcattcactcagtggaatctcagtagtaagccatgtggaatga