

ANGPT1-angiopoietin 1 Gene

Reference: PTXBC029406

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AGP1; AGPT; ANG1; angiopoietin-1; ANG-1; angiopoietin 1		

For research use only.

GENE DESCRIPTION	angiopoietin 1
GENE ID	284
ORIGIN SPECIES	Human
ACCESSIONS	BC029406
VECTOR	pDONR223
NCBI SYMBOL	ANGPT1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacacagtcacaacctgtcaatcttgcactaaagaagtgttactaaaggaggaaaaagagaggagaagaaccatttagagactgtgcagatgtatatcaagctggttttaataaaagtggaaatctacactatttatataataatgccagaacccaaaagggtgtttgcaatatgg
atgtcaatggggagggttgactgtaatacaacatcgtgaagatggagctagatttccaagaggctggaaaggaatataaaatgggttttggaaatccctccgggtgaatattggctgggaatgagttattttgccattaccagtcagaggcagtcacatgctaagaattgagtaattggact
gggaagggaaccgagcctattcacagtatgacagattccacatagggaatgaaaagcaaaactataggttaa