

MT1F-metallothionein 1F Gene

Reference: PTXBC029453

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MT1; metallothionein-1F; MT-1F; MT-1F; metallothionein 1F (functional); metallothionein-1F; metallothionein 1F		

For research use only.

GENE DESCRIPTION	metallothionein 1F
GENE ID	4494
ORIGIN SPECIES	Human
ACCESSIONS	BC029453
VECTOR	pDONR223
NCBI SYMBOL	MT1F
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacccaactgctcctgcgcgctgggtgtctctgcacctgcgctggttcctgcaagtgcgaagagtgcaaatgcacctcctgcaagaagagctgctgctctgctgccccgtgggctgtagcaagtgtgccagggtgtgtttgcaaagggcgctcagagaagtgacgtgctgcgactga