

MT1A-metallothionein 1A Gene

Reference: PTXBC029475

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MT1; MT1S; MTC; metallothionein-1A; MT-1A; MT-1A; metallothionein 1A (functional); metallothionein 1S; metallothionein-1A; metallothionein 1A		

For research use only.

GENE DESCRIPTION	metallothionein 1A
GENE ID	4489
ORIGIN SPECIES	Human
ACCESSIONS	BC029475
VECTOR	pDONR223
NCBI SYMBOL	MT1A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacccaactgctcctgcgcactgggtgctcctgcacctgcaatgcaagagtgcaaatgcaactcctgcaagaagagctgctcctgctgccccatgagctgtgccaagtgtgccagggtgcatctgcaaggggcatcagagaagtgcagctgctgtgcctga