

C9orf24-chromosome 9 open reading frame 24 Gene

Reference: PTXBC029484

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CBE1; NYD-SP22; SMRP1; bA573M23.4; spermatid-specific manchette-related protein 1; ciliated bronchial epithelial protein 1; testis development protein NYD-SP22; chromosome 9 open reading frame 24		

For research use only.

GENE DESCRIPTION	chromosome 9 open reading frame 24
GENE ID	84688
ORIGIN SPECIES	Human
ACCESSIONS	BC029484
VECTOR	pDONR223
NCBI SYMBOL	C9orf24
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttctgttctccgtaagaccaggaccctatcagtaacctacagtgactcctacagggtccacccctcaagagggtctataaggaccacccctgtgtgctgggaagccaacaagtttctgactcgggtctgactcataccatggagcgacatgtggatccgaagccctgcag
aagatggccaaatgtgtgtacaggactacacttataggggtccatatacaggccacccctacttgctgagaagtactggctttctcaaggaggagcagacaatgcagcccaactacctgggcagtgtgactgtacaacacatggaggatggaaccttaacacagcagctgtgcaac
aagtataccactaccttctcggctgcctaaaggaggccagatggagacagcagttcaggaatgccccttggaatgccctcctaggccggagcggctcaatgcctacgagcgcgaagtgtggtgaacatgctgaactcactgtcgcggaaccagcagctgccggtacacgcccc
gatcgggtgctgtgacccgtgcccggcgccctgcccctcatgtttacgaaagtcttgcggccgcccactactgtctgcgggatggaactactacgccagcggggcgccctgcaccgaccgcccctgcggccttgtgtccgggagcaacagactatgtgtacctccctacg
agcaccggcccgaatgcagtgctgtttacaactccccgcgtcactacccatccgaacctag