



Antibodies and Proteins

ProteoGenix SAS
Espace Européen de l'Entreprise
15 rue de La Haye
67300 Schiltigheim, France
(0)3 90 20 54 70 - Fax +33 (0)3 88
contact@proteogenix.fr
www.ProteoGenix.science

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88
contact@proteogenix.fr
www.ProteoGenix.science

TMEM120A-transmembrane protein 120A Gene

Reference: PTXBC029487

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NET29; TMPIT; transmembrane protein 120A; transmembrane protein induced by tumor necrosis factor alpha		

For research use only.

GENE DESCRIPTION	transmembrane protein 120A
GENE ID	83862
ORIGIN SPECIES	Human
ACCESSIONS	BC029487
VECTOR	pDONR223
NCBI SYMBOL	TMEM120A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcagccccgcccccggggccctctggggcactgcctgcgggacfcgggagatctacagcaggactctcagaacatccaggagaccatcggctctaccgcctgaagctcggaggagctgaccaaactcagaacaattgcacagctctcatcaccgggcagaagaagcgctcca
agctcctgcgcctcgccctgaagaataatgcaaaacctctcccaccagagcggccgacggggccgcacaggagctcgaagaacggctgaagaagcgcgcgaagcctctcttgcattgcagctgccttaattgccttaagaagaatgattgattctactgagctcgtctctgggaacgtcaacgctc
acgctcctgcagcagaacgaactaattgctctacagcagcagatgtagaattcgaactctctacatcatctctctctctcactctgcctctcctctcctcactcggctcctctcctcgaacccgggtgcacagatgctgcttaactctcctgctgactacacccctcaccctcggagagc
atctctatcaacaacggctcccggctggcagggcgccggcccgaggggggggagcatggaatgggggtgcgcagccctcgggaagcaggaacggcgctgggggtgatgggtgctcctgctcgcacgatcaaaaggctgggtgggtgctcatcactactgctccactctcctcggg
atgatcgtcagctggtgcgcagcgtctcatgatcaaaatctccggaaacaattctctctctctcatgatcagacgtcgtcagttctcagctactactacagagcgcctgcctctacccgctcggcgctggcgagcggcaacacatgggaactcactgtagggggtctcagtc
ctgtagtgcggggcctccactcctcgtccttctcttttttttggacactctgcgaactcttttaa

More info: Please contact us for any question: contact@proteogenix.fr