

## EYA3-eyes absent homolog 3 (Drosophila) Gene

Reference: PTXBC029500

|                |                                                                                         |                     |                |
|----------------|-----------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                     | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                    | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | eyes absent homolog 3; eyes absent 3; EYA transcriptional coactivator and phosphatase 3 |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | eyes absent homolog 3 (Drosophila)                                                                |
| <b>GENE ID</b>          | 2140                                                                                              |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC029500                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | EYA3                                                                                              |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atggaagaagagcaagatttaccagagcaaccagtgaagaaagccaagatgcaggaatcaggagagcaaaactataagtcagtaagcaatccagatgcagtgcagagaagcctgaacatcaagcctgcttcaaaccttccatgtcagaggaattatgacatgcaccgattacatcc  
ctcgtcatccaatgattatacctcacaatgtattctgcaaaacctatgcacatatctctcagttctgttttcggaactgcttaccctggacagactcaataccagacactacagcagactcaacctatgctgtctaccctcaggcaacccaacgtatggactaccttttgggtcattgtg  
gccaggatgaacacctgaaagtgttttaattcagactccatctcgaagtcacacagtggtttctacctgcactacaggggttaaccacaagccagccaagcccagcacattattcttacccttcaagcacaatgccagcctgatactacttcttacaattgccaattccagcag  
cagcagtagccagcatctcaaacaggattatccacctatactattcttggtcagaatcagtcaggcctgctacccagctccagctttggagtcacaggtcagactaacagtgatgcagagagcaccacattagcagcaaccacataccagtcggagaagcctagtgcatggcgct  
gcacctgcagcacagagacttctctgagacccttctacaagtcacattttgcccagactacacaaagtaagatactgatgcagtcaggaaaaacatgactagcaagaaccgggcaagaggaagctgagccacttctccagacagtgattagaacgggtatttctgtgg  
gacttggatgaacatcatcttccactacttcttactggtatcctatgccagaataatggaaggacccaacagtagtgattggctcagggttaacaatggaagaaatgattttgaagtggctgatactcatctattttcaatgacttagagacatag