

C6orf206-chromosome 6 open reading frame 206 Gene

Reference: PTXBC029519

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C6orf206; CILD12; MRPS18AL1; radial spoke head protein 9 homolog; radial spoke head 9 homolog		

For research use only.

GENE DESCRIPTION	chromosome 6 open reading frame 206
GENE ID	221421
ORIGIN SPECIES	Human
ACCESSIONS	BC029519
VECTOR	pDONR223
NCBI SYMBOL	C6orf206
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacgccgacagcctcctgctgtctctggagctggcgtccggcagtgaggcagggcctcagcccgacgctcggcctcgtctcagctctcttatgctgttaagcgcgactaccgctatgatcgggttctctctggggccgcatccttggcctcgtcgcgattactacatcgcgca
gggcctgagtgaggaccagctcgcaccgcgcaagacgtctatagcctgaactgcacagagtgaggcctctgccccctgccacagaggagatggtggcgagctcgtctgtggaaggccgcttcatggggacccatatacgaatgaacacactgagctgcagaagtgat
gaaggtgaaaaagtcttgaagaagaatagtggtccagatcaaggagagaccgcttgggtgtctgctgaccagattgacaaggctgtggccatcatccccgagggcgcctcttcaagacccctttggacccacccatgtcaatcggaccttgaaggactgtcctgtctgaggcc
aagaagctcagctcctacttccattcaggagcctgttgagctaaagaagaacgttggagaaggctgacgtggacccctcctggtattcatgactccttgagcatgacattccaaagggtcctggagcatccagatggagaggggcaatgccctggtgtgtcgcgcagcct
gctctggcgggcctcacttaccatgctccccgcaccaagaactatggctacgtctacgtgggcactggcgagaagaacatggactggccttcatgctatag