



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

MGC33894-transcript expressed during hematopoiesis 2 Gene

Reference: PTXBC029527

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C17orf103; Gtlf3b; protein NATD1; N-acetyltransferase domain-containing protein 1; gene trap locus F3b; protein GTLTF3B; transcript expressed during hematopoiesis 2; N-acetyltransferase domain containing 1		

For research use only.

GENE DESCRIPTION	transcript expressed during hematopoiesis 2
GENE ID	256302
ORIGIN SPECIES	Human
ACCESSIONS	BC029527
VECTOR	pDONR223
NCBI SYMBOL	MGC33894
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgcactcggcgccgcgctggcgcgtggcgcgtggagcaggcgtgcccatcgcgtggagcagcagccggcgccgcagctactgtccggctcaacggtgtcatgaccgggctctctgtcatgtagctacgtgggcaagcggalctgfgacctgcagcacaccga
gtgccagatgctacagctggcgctggcgtccaaagcactctgcaagcggccgccgtgactctgtgtggaggaagacttgaagcccatctaccctctgtatcatcagaagtactctgaagaaacccccctgccagttacctgtgaagcgcctgcagccgttaa