



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

C19orf41-chromosome 19 open reading frame 41 Gene

Reference: PTXBC029535

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C19orf41; PLAL6978; PRO21961; SCRL; izumo sperm-egg fusion protein 2; IZUMO family member 2		

For research use only.

GENE DESCRIPTION	chromosome 19 open reading frame 41
GENE ID	126123
ORIGIN SPECIES	Human
ACCESSIONS	BC029535
VECTOR	pDONR223
NCBI SYMBOL	C19orf41
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]