

GLIPR1L2-GLI pathogenesis-related 1 like 2 Gene

Reference: PTXBC029557

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GLIPR1-like protein 2; GLI pathogenesis related 1 like 2		
For research use only.			

GENE DESCRIPTION	GLI pathogenesis-related 1 like 2
GENE ID	144321
ORIGIN SPECIES	Human
ACCESSIONS	BC029557
VECTOR	pDONR223
NCBI SYMBOL	GLIPR1L2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaggccgcaaggccctcgcggggagtgaggggccagtcctacccctggcagtagggggcgtttgaagctcggtctgtgagctgtggctactgctactgggttctagttgaacgccagattttgccagacagaggagacgtagactttatcaacgagtagctgaacctcc
acaatgagctcgggggcagctcattccccgagggtctaacttgcgcttcagcttgggatgtagctttatcacgactgctagagcatggggaaaaaatgtttgttacgcataatattttacaagatgtacaaatggccatccataattttatgtattggtgaaaaatgtgggtcggcc
ctgaaaatgaattactgcaagtattgctatcagaagttggcatgcagagaagaaaaatgtacaattttgaaaatggcagttgctctggagactgttctaatattacagctgtttgggaccactcttacaagttggtgtgctgtactccatgttcaaaaatggacatattatcatgcagcaatttt
catatgcaactatgcggcaggaggaactgacgagaagacattatgaaccaggaattttgtactcgtatggcagacgtgacaaatgcacagattttctatgcagtaagataaagaaaataacatgaaaaaatgcataatggattggacaagaaaaataagcgattgaacactagttt
ttatggctcatgtaa