

ZFAND2A-zinc finger, AN1-type domain 2A Gene

Reference: PTXBC029558

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AIRAP; AN1-type zinc finger protein 2A; arsenite inducible RNA associated protein; zinc finger, AN1-type domain 2A; zinc finger AN1-type containing 2A		

For research use only.

GENE DESCRIPTION	zinc finger, AN1-type domain 2A
GENE ID	90637
ORIGIN SPECIES	Human
ACCESSIONS	BC029558
VECTOR	pDONR223
NCBI SYMBOL	ZFAND2A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagtttctgatttgggaagcattgttcagaaaagactgcaagcagctagattttctcagtaaaatgtgcatgtaaacagattctgtaaaagatcttccatagctgcacataagtgctcgtttgcatccagaaggatgtgctgctgccagatgccactctgtaatacccaa
tcccagtaaaaaggccagataccagacgtggtgttgatcacattgacagagactgtgactctcaccctgggaagaagaagaagaattttacataccgttgctcaaaagaggctgcaagaagaagaagatgctgcagatgggtatgtgccaatgtcacggcaactctgtatc
cagcacagacaccttggaccacagctgcagacacgggagtcgccccacatcaagctgggtga