

VPS53-vacuolar protein sorting 53 homolog (S. cerevisiae) Gene

Reference: PTXBC029560

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	VPS53, GARP complex subunit; HCCS1; PCH2E; hVps53L; pp13624; vacuolar protein sorting-associated protein 53 homolog		

For research use only.

GENE DESCRIPTION	vacuolar protein sorting 53 homolog (S. cerevisiae)
GENE ID	55275
ORIGIN SPECIES	Human
ACCESSIONS	BC029560
VECTOR	pDONR223
NCBI SYMBOL	VPS53
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagtttctgtctggacacccactcgtgaagatggtcctctcgtcgtccctccatcagctcgcaggtggtgaggaaggcaccgccagctacaccaagatcgtgtcaaaaggcatgccgggctgagatgatcctcaaggtagtgatggccctcatgaaccgttggtggtt
tgftgacaactacatcaaaattctcacagactgcaacacagaaacctttcagaagatactggacatgaagggtctgaagaggagtgagcagagcagcatgctggaactcctgcgcagcggctccccgcaccgccctcggggcagaaagctcggctcactgtccctgacggcgccg
acaccagagcaagagtcgtcacgcatccgcaagctcgagaaactcaaaaaagagactgtag