

C8orf56-chromosome 8 open reading frame 56 Gene

Reference: PTXBC029562

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	chromosome 8 open reading frame 56		
For research use only.			

GENE DESCRIPTION	chromosome 8 open reading frame 56
GENE ID	157556
ORIGIN SPECIES	Human
ACCESSIONS	BC029562
VECTOR	pDONR223
NCBI SYMBOL	C8orf56
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtccttaagtctctggcatcctcagagcaagacaaaaagatgggggcaagcgaaggttaacccacgtggggcctctggaagtatggagggccctctcttctctcttccacctctggccagcgagggcagaactcacaggcaaacctgggttctgcacagatgcagctgcatttt
gaacctagaggaaagcatggacagcgactgggggcttgggtgggggtgtcgtccctagaaggcaccctttctgatatgggtcagatgggccttgggtgcacacaggcagcgtccaggatggggacattga