



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

MGC39584-hypothetical gene supported by BC029568 Gene

Reference: PTXBC029568

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

uncharacterized LOC441058

For research use only.

GENE DESCRIPTION	hypothetical gene supported by BC029568
GENE ID	441058
ORIGIN SPECIES	Human
ACCESSIONS	BC029568
VECTOR	pDONR223
NCBI SYMBOL	MGC39584
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatgaaacccaggcggagacgggggaagcagcacgggatcccagcctcaggcctgcacggacggtgttggttggtcacttcatggctaagaaagcaactccaaccacacactgcaaataatgaaggaccacgaaaggcacgcaggctctggttgcatag