

LOC100130950-similar to hCG1991536 Gene

Reference: PTXBC029580

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized LOC100130950		
For research use only.			

GENE DESCRIPTION	similar to hCG1991536
GENE ID	100130950
ORIGIN SPECIES	Human
ACCESSIONS	BC029580
VECTOR	pDONR223
NCBI SYMBOL	LOC100130950
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcaccaaggaaaggccatttgggggacatagcaagaaggcagctgtgtgcaagccaggaggagaacctcaccagaaacctcacagaaatcgcccaacacctcatctcgggcttcagctcagaacaacgagaaaacaaattctgtgtgtcctggttccatttgaattggca
ccctaataaaccgttttctgttagccctggaaatggctgaaatgaggggaggggaataaagtgtctataattcagagcatggagcacaccttctgttcacctgagctggattcggagaggttaa