

ADIG-adipogenin Gene

Reference: PTXBC029594

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SMAF1; adipogenin; adipogenesis associated; small adipocyte factor 1 (SMAF1)		
For research use only.			

GENE DESCRIPTION	adipogenin
GENE ID	149685
ORIGIN SPECIES	Human
ACCESSIONS	BC029594
VECTOR	pDONR223
NCBI SYMBOL	ADIG
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgactccagtgtgtgcttggattggagccctggagcaaaagccagctgagtttctggaaggggacactccacggccaagagaaggagaggcctgctgtgagcctgtgtgccaggtgaggcactccaggggagggggagcctcaaggcccaacaaagccttgggc
agctgctatgtgggcaagaggctgcctccaccatattagatttcttctgagtcagcagggcagtggaatggcagaagtgatgggagagacttgcaggaggcaagaggactttggcaactgtgtgtgcctccctggaaccccaactcctcctctcagaccgacatgtga