

IQCF1-IQ motif containing F1 Gene

Reference: PTXBC029595

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	IQ domain-containing protein F1; IQ motif containing F1		
For research use only.			

GENE DESCRIPTION	IQ motif containing F1
GENE ID	132141
ORIGIN SPECIES	Human
ACCESSIONS	BC029595
VECTOR	pDONR223
NCBI SYMBOL	IQCF1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaggagaagcagccccaaagacgaaggaaacctcaaaagaagatgagctcagcagagaaggatgccaacctttgtccttaggagcagagtcgaaaggcagaggctaaactccagttctgttgagacacagacagtggaatgccaatgagaaatcagaaaaaccccc
agaaaacaaaagaagctctctgacaagatacggtagccaccaagatcaaggcctgtgtggcggggcaccctggtgctgctgggcactgttgacgcagccctcagtgctgcatcattcagtgctgtggcgctgatactgtccaagattctgaagaaggcggcaggcagcgcta
gaggccttctccgggaaggagtgaggcagcagtcacactgcagtcagccaggcccgcatgtggcgcatccgcagacgctattgccaggtgtcgaatgtctgcatcatccaggcttactggaggtgccgctcctgtgtctccgggggttcataaggccagtcacagtcacagccaac
cagctgcatctcagctggagatctgtgtgactcaggcccttcattgtgacagagtgattcccttcaataaaggaaatga