

SYT9-synaptotagmin IX Gene

Reference: PTXBC029605

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	synaptotagmin-9; synaptotagmin IX; sytIX; synaptotagmin 9		
For research use only.			

GENE DESCRIPTION	synaptotagmin IX
GENE ID	143425
ORIGIN SPECIES	Human
ACCESSIONS	BC029605
VECTOR	pDONR223
NCBI SYMBOL	SYT9
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctggcgaaggtggtggaggagatttagctttcaaggagcagagatctcagtgagcctgtgacccttgggtgactgctgtgtctcgtctctttggcgtgtctctctctgtatcttgaaactctgtgggtccgtggcgagaacgagccctgcccctggttagcaagacaa
caaccaggagcccccttaactacatggacacagagaccaatgagcaggagaacagtgaggacttctagacctcccacgcccctgcccctgactcctccatgaagatcagccacacctcccctgacattccccctccaccagacgggatccaggaactgtgcccgtgctccgc
gtgcagcgcgaagtcacagagccaaactcgtgcccgcataattcaatccgaagacaactcaacttgcacaaacccggacttcaatatccagcagctcaaaaacagggaacagttgactggaattggtagaattaaaccagagttatataagcagaggtcattggataatgatgacggga
gacggagttaacagcaaggctgtgtgggaaactgaacttcattttaaaatagactgtgacttagcagctcagtagtagaagattcacaaagctgtcaatttggccccaaggacttttctgggacttcagatccttatgcaagatctatttctcctgacgggaaacaaacaccagactaaagt
tcacagaaagaccctgaacctgtgtttgatgaagtgtttatttccgttccctacaatgaccttgaagcacggaagcttcaactctctgtgtacgacttgcagaggttctctcgtcatgacttaatcggccaagtgggtggatcacttctagacttggctgatttccccaggagtgatccttt
ggaaggatcgaatatgtcaccaatatcttagtag