

LOC554206-hypothetical LOC554206 Gene

Reference: PTXBC029609

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical LOC554206
GENE ID	554206
ORIGIN SPECIES	Human
ACCESSIONS	BC029609
VECTOR	pDONR223
NCBI SYMBOL	LOC554206
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctgctggacacagctccccgaacctccacgcccattggccactagacagaggagctctcttccacctggaattgcgacgctggcgacgagggcatggctgcacctgcaagatgcttctgtgcaagaagcgcctgtcgggcgctggattcggggct
ggcatctgggacgctggctgaggtgggagcgggctgcatctgaagaatacgtcggaggtggcaggttgcctccccgcctcgcacgacctcgttccacctgtgaaatgcacagaacagggttcatttattaacgaatcgtttctgagctcctgtgtgagccagcttggag
caagcctgggtactgtggtgggagcagcccctga