



Reference: PTXBC029653

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ATBF1; ATBT; ZNF927; zinc finger homeobox protein 3; AT motif-binding factor 1; AT-binding transcription factor 1; ZFH-3; alpha-fetoprotein enhancer binding protein; zinc finger homeodomain protein 3; zinc finger homeobox 3		

For research use only.

GENE DESCRIPTION	zinc finger homeobox 3
GENE ID	463
ORIGIN SPECIES	Human
ACCESSIONS	BC029653
VECTOR	pDONR223
NCBI SYMBOL	ZFHX3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]