

CST9L-cystatin 9-like Gene

Reference: PTXBC029656

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CTES7B; bA218C14.1; cystatin-9-like; testatin; testicular tissue protein Li 45; cystatin 9-like		

For research use only.

GENE DESCRIPTION	cystatin 9-like
GENE ID	128821
ORIGIN SPECIES	Human
ACCESSIONS	BC029656
VECTOR	pDONR223
NCBI SYMBOL	CST9L
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgggctgcccgtggaaggaggtctgtcctggcgctgctgtcttcttaggtcccagatcctgctgatctatgcctggcattccacagcaaggagctgtgatgaacacaatgtcatggctcgttacctccctgccacagtggagtttctgtccacattcaaccaacaga
gcaaggactactatgcctacagactggggcacatctgaattcctggaaggagcaggtggagtccaagactgtattctcaatggagctactgctggggagaactagggtgtggaaatttgaagacacattgacaactgccattccaagaagcacagagctgaacaatacttccactgt
tcttcaccatcagcaccaggccctggatgactcagttcagcctcctgaacaagacctgctggagggaattccactga