

TMEM61-transmembrane protein 61 Gene

Reference: PTXBC029775

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	transmembrane protein 61		
For research use only.			

GENE DESCRIPTION	transmembrane protein 61
GENE ID	199964
ORIGIN SPECIES	Human
ACCESSIONS	BC029775
VECTOR	pDONR223
NCBI SYMBOL	TMEM61
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccctgccccagatgtgtgacgggagccacttggcctccacccctcgtattgcatgacagtcagcggcagcagtggtttctgtggcgggacgctctgttctgttggaggcgaaggggatgcaaccgcccagctggccagctggccccacccacggagtatccggtgcctg
agggccccagccccctgtcctcaggtccgtcagcttctgtctgtcgttgcaggtggcctgtctgtctcattggcctgtgtgttccgtcaaggccagatcccagggccacctcgtgggacccctatcacctctccagagacctgtactactcactgtggagtccagagaaggagag
ctgcaggacccccaaagtgggtgacatccccacttacgagggaagccgtgagcttcccagtgccgagggggcccccaacaccacctgcataccctacggagggaagccctggagccaagtggatcaggggatgccctgtcagcaccagcccgcctggcctcccccagctatgaga
gcatcagcctgtcttgcgttctgtcagagacgacaccgagtgccacacgctcctgtcaggcctgttgcagactgcacgggaggaagtaa