

DMRTC1-DMRT-like family C1 Gene

Reference: PTXBC029799

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	doublesex- and mab-3-related transcription factor C1; doublesex-mab-3; DMRT like family C1		

For research use only.

GENE DESCRIPTION	DMRT-like family C1
GENE ID	63947
ORIGIN SPECIES	Human
ACCESSIONS	BC029799
VECTOR	pDONR223
NCBI SYMBOL	DMRTC1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctgccccccaaagctccatccgtgtcaggaatttgaccatcagagcaggagccctcactgggaaggagaaacaacatgtgcagcccgagacccacatcttcacagccccgaggaggggagctccaaggggctctgctgtggccaggccccagaacctttgtctgcc
ctgtactccagtgacctggagcagcaactggtttctcttctgggatccccacagggccctgcctgcccagcatatgtcacaactctgatctccagccctgtgccaccttgacctcttctactgcagccacagggtctgggaaagtag