

ADORA3-adenosine A3 receptor Gene

Reference: PTXBC029831

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	A3AR; adenosine receptor A3; adenosine A3 receptor		
For research use only.			

GENE DESCRIPTION	adenosine A3 receptor
GENE ID	140
ORIGIN SPECIES	Human
ACCESSIONS	BC029831
VECTOR	pDONR223
NCBI SYMBOL	ADORA3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcccaacaacagcactgctctgcatggccaatgttacatcacaccatgaaattttcattggactctgcccatagtgggcaacgtgctgcatctgctggaagctgaacccagcctgcagaccaccacctctatttcattgtctcttagccctggctgacattgctgttgggt
gctggtcatgcttggccattgtgtcagcctgggcacacaaatccattctacagctgctttttatgactgctactgcttatctttaccacgctccatcatgctctgctggccatgctgtggaccgatacttgcgggtcaagcttaccgtcagatacaagagggtcaccactcacagaa
gaatatggctggccctggccttggctggtgtcattcctggtggatggaccccatgttggctggaacatgaaactgacctcagatgaccacagaaatgtcacttcttcatgccaattgttccgtcatgagaatggactacatggtatctcagcttccctcactggatttcatccc
cctggtgtcatgtgcacatctatctgacatctttacatcattcggaaacaaactcagctgaacttcttaactcctaaagagacaggtgcatcttattggacgggagttcaagacggctaaagtcctgtttctggttcttctgttctgtctgtcatggctgctttatctatcatcaactgcatcatcta
cttaatggtgaggtaccacagcttgtgtgtacatgggcacctgctgtccatgccaactcatgatgaacctatgctctatgctataaaataaagaagtcaaggaaacctacctctgacatcaaaagcctgtgtgtgtgctgccatccctctgattcttggacacaagcattgagaagaatt
ctgagtag