



Reference: PTXBC029854

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	INT1L1; IRP; protein Wnt-2; Int-1-related protein; int-1-like protein 1; secreted growth factor; wingless-type MMTV integration site family member 2; Wnt family member 2		

For research use only.

GENE DESCRIPTION	wingless-type MMTV integration site family member 2
GENE ID	7472
ORIGIN SPECIES	Human
ACCESSIONS	BC029854
VECTOR	pDONR223
NCBI SYMBOL	WNT2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]