



Antibodies and Proteins

ProteoGenix SAS
Espace Européen de l'Entreprise
15 rue de La Haye
67300 Schiltigheim, France
(0)3 90 20 54 70 - Fax +33 (0)3 88
contact@proteogenix.fr
www.ProteoGenix.science

ENTPD3-ectonucleoside triphosphate diphosphohydrolase 3 Gene

Reference: PTXBC029869

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD39L3; HB6; NTPDase-3; ectonucleoside triphosphate diphosphohydrolase 3; CD39 antigen-like 3; NTPDase 3; ecto-ATP diphosphohydrolase 3; ecto-ATPDase 3; ecto-ATPase 3; ecto-apyrase 3		

For research use only.

GENE DESCRIPTION	ectonucleoside triphosphate diphosphohydrolase 3
GENE ID	956
ORIGIN SPECIES	Human
ACCESSIONS	BC029869
VECTOR	pDONR223
NCBI SYMBOL	ENTPD3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr