



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

C13orf16-chromosome 13 open reading frame 16 Gene

Reference: PTXBC029889

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C13orf16; bA474D23.1; testis-expressed sequence 29 protein; testis expressed 29		
For research use only.			

GENE DESCRIPTION	chromosome 13 open reading frame 16
GENE ID	121793
ORIGIN SPECIES	Human
ACCESSIONS	BC029889
VECTOR	pDONR223
NCBI SYMBOL	C13orf16
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaataacgtctggaagtgaaagaacctccgcgcacacctctgaagcaattcacagttgtgtgatgttccctctgatgaacatttgtaactacaacgtctccaggagcagatgccaggagctcgggtgctgcttctacgaagcgctgtctacaagaagcgggttccattatcaccaagtggtctctcgtgattgatcatcctgctggcgccctgtcatcacatcatctacagagattctcaggagagacaggaaagaaagccatccctctgtgatgtgctctccagcagatgagctcagcgaaaagcggcgctctctacagcagcaatgttagggctgaagcctgcgagtctcggccctcagatgaagagtgatgacgaggatgaagatgatgaacgacaaataacagaagccgaagaaacatgagactga