

TFEC-transcription factor EC Gene

Reference: PTXBC029891

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TFEC-L; TCFEC; TFE-C; TFECL; bHLHe34; hTFEC-L; transcription factor EC; class E basic helix-loop-helix protein 34		

For research use only.

GENE DESCRIPTION	transcription factor EC
GENE ID	22797
ORIGIN SPECIES	Human
ACCESSIONS	BC029891
VECTOR	pDONR223
NCBI SYMBOL	TFEC
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacccttgatcatcatcatcaatccttaaatgtgcacaacctgcagtgccaagtggggcctctgtgcagcatgcacacacaactctggacagtgtgctgcctcacagaaaaccctcaccaagttactagctattgggaagaagatgacaatgcacaatggcatttat
ctggagatattttggtgtatagcggfgaacaagggaatttcaccaattaacatgggcttacaagtgtcttctgtccaagtgtctaccaatgaaaagagaaattacagaaactgacactagagcttttagcaaaagagagacaaaaaaggacaaccacaacctcattgaaagaagaaga
gggtataatattacccaatcaaggagcttggcactcttattccaaagtctaagatcctgatatgcgctggacaagaaggaaccttctaaagcatcagtgagtagcatcaagtggtctacaaaaagaacaacagagagcccgagaattggaacacagacagaagaattagagcaggct
aacaggcgacttctacttgcgattcaggttttataagaatgtag