

TSPAN16-tetraspanin 16 Gene

Reference: PTXBC029908

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TM-8; TM4-B; TM4SF16; tetraspanin-16; tetraspanin TM4-B; transmembrane 4 superfamily member 16; tspan-16; tetraspanin 16		

For research use only.

GENE DESCRIPTION	tetraspanin 16
GENE ID	26526
ORIGIN SPECIES	Human
ACCESSIONS	BC029908
VECTOR	pDONR223
NCBI SYMBOL	TSPAN16
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctgaatccacactccgtattcttcctgaagaactgttatctttactcaatggctcgtgctgtctgcatcactagttggcctggcattggtgtaaatgtggaggggcctctgacgaatgtctcggcgtgctccgcatacctcctcagttggcaacctgtgcctggtg
atgggatgcatcacgggtactccttggctgtccgggtggtatggagcgactaaagagagcagaggcacgctctgttttgcacctcgtcaatggttattgtctcactatggaagttacagctgccacagtggctcctcttcttccaattgttgagatgtgaccttggaacacaccttcgtgac
cctgaggaagaattacagaggttacacagagccagacgactattctacacagtgggaacttggtcatggagaaggactcctccaagtatttctctctctctatag